

Gastrointestinal Panel

QAPs™ for Molecular (MDx) and Immunoassay (EIA) Diagnostics



- Formulated with the pathogen's whole genome (MDx) and/or the target antigen (EIA)
- Whole workflow controls (including elution and/or extraction)
- Samples meet clinical laboratory biosafety requirements
- Liquid formulations are stable for 1+ years when stored at 8° C; 1mL and 0.5mL volumes
- Swab formulations are stable for 1+ years when stored at 2-30°C
- Commutable and cross-platform compatible
- Verified performance on leading IVD platforms

Compatible With

Product	Format	PTDx™	PROCEEDx™	REDx™	MDx	Immuno
Single Analyte Products						
<i>Cryptosporidium parvum</i> Positive	Vial	PT-37-01	VP-37-01	RED-37-01**	✓*	✓
<i>Giardia lamblia</i> Positive	Vial	PT-38-01	VP-38-01	RED-38-01**	✓*	✓
Rotavirus Positive	Vial	PT-35-01	VP-35-01	RED-35-01**	✓*	✓
<i>C.difficile</i> GDH Positive Swab	Swab	PT-S-39-01	VP-S-39-01	Coming soon	x	✓
<i>C.difficile</i> GDH/Toxin A/ Toxin B Positive Swab	Swab	PT-S-39-02	VP-S-39-02	Coming soon	x	✓
Campylobacter OMP Positive Swab	Swab	PT-S-40-01	VP-S-40-01	Coming soon	x	✓
Multiple Analyte Products						
<i>Cryptosporidium parvum</i> + <i>Giardia lamblia</i> Mixed Positive	Vial	PT-37-02	VP-37-02	RED-37-02**	✓*	✓
Negative Analyte Products						
Negative Sample	Vial	PT-99-01	VP-99-01	Coming soon	✓	✓

* Products are QC- released with immunoassays but are compatible with molecular assays

** Products are only for sale in Canada and the United States



Intended Use

- REDx™ products are unassayed controls that can be used for ongoing External Quality Control (QC), benchmarking laboratory testing performance, operator training & competency and end-to-end workflow evaluation.
- PROCEEDx™ products are used for assay verification/validation, training, and R&D purposes. Products are for RESEARCH USE ONLY!
- PTDx™ products are used as challenge samples in laboratory proficiency testing (PT) and external quality assessment (EQA) programs

Need to Know More or Want to Order?

Web: microbix.com **Email:** customer.service@microbix.com

Phone: 800-794-6694 toll free in North America or +1-905-361-8910 worldwide



Microbix Biosystems Inc.
265 Watline Ave. Mississauga,
ON Canada L4Z 1P3
+1-905-361-8910
www.microbix.com

www.microbix.com