

Genital Ulcer Disease Panel

QAPs™ for Molecular and Immunoassay Diagnostics



- Formulated with the target pathogen's whole genome
- Whole workflow controls (including extraction)
- Inactivated to meet clinical laboratory biosafety requirements
- Commutable and cross-platform compatible
- Liquid formulations are stable for 24 weeks when stored at 2-8°C; 1mL volume
- Copan FLOQSwab® formulations are stable for 1-2 years when stored at 2-30°C
- Verified performance on leading IVD platforms

Compatible With

Product	Format	PTDx™	PROCEEDx™	REDx™	MDx	Immuno
Single Analyte Products						
HSV 1 Positive	Vial	PT-02-M1	VP-02-M1	RED-02-M1*	✓	✓**
HSV 2 Positive	Vial	PT-23-M1	VP-23-M1	RED-23-M1*	✓	✓**
HSV 1 Positive Swab	Swab	PT-S-02-01	VP-S-02-01	Coming soon	✓	✓**
HSV 2 Positive Swab	Swab	PT-S-23-01	VP-S-23-01	Coming soon	✓	✓**
<i>T. pallidum</i> Positive Swab	Swab	PT-S-80-01	VP-S-80-01	Coming soon	✓	X
Varicella Zoster Virus Positive Swab	Swab	PT-S-03-01	VP-S-03-01	Coming soon	✓	✓**
Monkeypox West African Positive Swab	Swab	PT-S-81-02	VP-S-81-02	Coming soon	✓	X
Multiple Analyte Products						
HSV1&2/VZV/Syphilis Positive Swab	Swab	PT-S-03-M1	VP-S-03-M1	RED-S-03-M1	✓	X
HSV 1 & 2 Positive	Vial	PT-02-M2	VP-02-M2	Coming soon	✓	✓**
HSV 1 & 2 Positive Swab	Swab	PT-S-02-M2	VP-S-02-M2	Coming soon	✓	✓**

Note: * Products are only for sales in Canada and United States

Note: ** Products are QC- released with molecular assays but are compatible with immunoassays

Intended Use



- REDx™ products are unassayed controls that can be used for ongoing External Quality Control (QC), benchmarking laboratory testing performance, operator training & competency and end-to-end workflow evaluation.
- PROCEEDx™ products are used for assay verification/validation, training, and R&D purposes. Products are for RESEARCH USE ONLY!
- PTDx™ products are used as challenge samples in laboratory proficiency testing (PT) and external quality assessment (EQA) programs

Need to Know More or Want to Order?

Web: microbix.com **Email:** customer.service@microbix.com

Phone: 800-794-6694 toll free in North America or +1-905-361-8910 worldwide



Microbix Biosystems Inc.
265 Watline Ave. Mississauga,
ON Canada L4Z 1P3
+1-905-361-8910
www.microbix.com

www.microbix.com